homecraft



STANDARD UTILITY CHAIR

Product Code: ZHI-ZH026

Adjustable height, contoured backrest and padded arm rests, vinyl upholstery and powder coated steel frame.

Overall Width	580mm
Seat Width	455mm
Seat Depth	460mm
Seat Height	450 - 605mm
Back Height	465mm
Max. User Weight	160kg
Warranty	12 months
Colour	Teal/Blue, Mocha

Complies with AS/NZS 4688.2:2000 or 4688.3:2000.



HIGH BACK CHAIR

Product Code: ZHI-ZH032

High back chair, adjustable height, contoured backrest and padded arm rests, vinyl upholstery and powder coated steel frame.

Overall Width	640mm
Seat Width	490mm
Seat Depth	480mm
Seat Height	455 - 605mm
Back Height	730mm
Max. User Weight	160kg
Warranty	12 months
Colour	Teal/Blue, Mocha

Complies with AS/NZS 4688.2:2000 or 4688.3:2000.



LEGREST STOOL

Product Code: ZHI-ZHSTOOL

The adjustable height legrest stool is designed with a 10° angle for comfort and support and should be placed with the higher end against the chair. The legrest stool matches the high back chair (ZHI-ZH032).

Seat Pad	400mm x 500mm x 40mm
Height Adjustable	420 - 570mm
Max. User Weight	160kg
Warranty	12 months
Colour	Teal/Blue, Mocha

ASSEMBLY

1. The extension legs are placed inside the tube and adjusted to height by depressing the snap button until it locks into the desired height hole.
2. Ensure that the snap buttons fully protrude through the same respective adjustment hole of each leg extension. This ensures that the leg extension are securely locked in position and an even height adjustment is achieved.
3. After adjustment and BEFORE every use, ensure that the seating platform is level and stable on a flat surface before using.
4. For HIGH BACK CHAIR only: After assembly of the leg base, remove the pin clips on the side tubes of the base, insert the backrest panel, line up the holes and insert the pin clips to secure the backrest in position.

Note: A pair of Extension Legs with Castors can be purchased separately to replace the extension legs at the back of the chair, making it easier to move the chair. To move the chair, tilt the chair back until the castors touch the ground - **DO NOT** tilt the chair or move the chair while someone is sitting in the chair.

CARE INSTRUCTIONS

Clean with disenfectant or mild detergent and water. Wipe chair to remove any soap residue. Use of caustic or abrasive cleaners will void warranty.

Manufactured for

Performance Health ANZ
Unit 3 / 3 Basalt Road
Pemulwuy NSW 2145
AUSTRALIA

Tel: 1300 473 422
Fax: 1300 766 473
sales@performhealth.com.au
www.performhealth.com.au

PERFORMANCE HEALTH™

Made in China