


2-in-1 Rollator and Transit Chair

Comfort when you need it

The 2-in-1 Rollator and Transit Chair offers the comfort you need for various uses. Each rollator includes a padded seat and backrest, an under-seat bag for storage and flip-up footrests to provide extra support.

Key features

- Suitable for use as a rollator or transit chair
- Lockable loop brakes for security and convenience
- Folds easily for transport and storage
- ISO 11199-2:2005 / Rollator Standard



A	Handle height 80 - 93cm (31½ - 36½")
B	Pushing handle height 90.5cm (35½")
C	Overall depth 74cm (29¼")
D	Overall width 58.5cm (23")
E	Seat height 56cm (22")
F	Wheel diameter 19cm (7.5")
	Seat width 35cm (13¾")
	Seat depth 33cm (13")
	Product weight 8.8kg (19.4 lbs)

Part#	Description
DAYS-2IN1-RED	Red
DAYS-2IN1-SIL	Silver - Gun Metal



Available colours

Red

Silver
Gun Metal

160 kg

25 st 352¾ lbs

Maximum user weight

1 YEAR

Product warranty

2 YEAR

Main frame warranty