



Alpha 426 Rollator

Features:

- lightweight aluminium frame
- unique height adjustable seat
- comfortable & durable PU seat
- adjustable handle & loop brakes
- 6" wheels & nylon storage bag
- easy to fold for travel or storage
- user weight capacity of 125kg

5 colours



ALPHA 426 ROLLATOR	Metric	Imperial
Approximate Total Product / Unit Weight	6.0 kg	13.2 lbs
Overall Width (product or frame when opened)	610 mm	24.0 inch
Overall Depth (product or frame when opened)	740 mm	29.1 inch
Width / Depth (product or frame when closed)	130 mm	5.1 inch
Overall Height (or height range if adjustable)	730 - 980 mm	28.7 - 38.6 inch
Armrest Width (outside dimension) Armrest	N/A	N/A
Width (inside dimension)	440 mm	17.3 inch
Armrest Height (from top of seat)	N/A	N/A
Back / Back-Rest Height (from top of seat)	270 mm	10.6 inch
Seat Height (from top of seat to floor)	460 - 580 mm	18.1 - 22.8 inch
Seat Width (actual seat, not available space)	350 mm	13.8 inch
Seat Depth (actual seat, not available space)	300 mm	11.8 inch
Wheels (dimensions of front)	152 mm	6.0 inch
Wheels (dimensions of rear) Recommended	152 mm	6.0 inch
Maximum User Weight	125 kg	275.6 lbs

