

## Alpha 427 Rollator

### Features:

- lightweight aluminium frame
- multi-height adjustable PU seat
- fold-up & removable backrest
- adjustable handle & loop brakes
- 7" wheels & nylon storage bag
- easy to fold for travel or storage
- user weight capacity of 125kg



Alpha 427 Rollator	Metric	Imperial
Approximate Total Product / Unit Weight	6.5 kg	14.3 lbs
Overall Width (product or frame when opened)	610 mm	24.0 inch
Overall Depth (product or frame when opened)	750 mm	29.5 inch
Width / Depth (product or frame when closed)	160 mm	6.3 inch
Overall Height (or height range if adjustable)	750 - 990 mm	29.5 - 39.0 inch
Armrest Width (outside dimension)	N / A	N / A
Armrest Width (inside dimension)	440 mm	17.3 inch
Armrest Height (from top of seat)	N / A	N / A
Back / Back-Rest Height (from top of seat)	270 mm	10.6 inch
Seat Height (from top of seat to floor)	460 - 590 mm	18.1 - 23.2 inch
Seat Width (actual seat, not available space)	350 mm	13.8 inch
Seat Depth (actual seat, not available space)	300 mm	11.8 inch
Wheels (dimensions of front)	178 mm	7.0 inch
Wheels (dimensions of rear)	178 mm	7.0 inch
Recommended Maximum User Weight	125 kg	275.6 lbs

