



The  
**Mobility Store**  
Solutions for Active Living

# Newsletter

**August 2025**

## Hello from The Mobility Store!

In the coming weeks, we're pleased to be expanding our Icare range with three new products - two innovative lift chairs and a newly designed bed. The new lift chairs are ideal for clients with complex or evolving needs, offering enhanced support and adaptability. The new IC333 Curvex bed features curved edges to promote greater safety and comfort for your clients. *Stay tuned for further details!*

### Icare VMotion Rise and Recline Lift Chair

The VMotion chairs fits into any home or clinical facility. It blends advanced functionality with everyday comfort.

#### Features

- Four Quiet Drive Actuators
- Pressure Care Inseat Options - Soft, Spring, Firm
- Ultra large dual storage side pocket
- Lumbar support can be rotated and vertically supported
- Available in fabric and vinyl



"With three interchangeable inset options, this bed offers flexible support to accommodate your client's evolving pressure care needs over time."

*Sue Beardsell, TMS*

#### Specifications

- 170kg weight limit
- Seat size - 500mm wide x 500mm deep
- Hand piece - back lit. Wired stretch cord
- Seat height - 500mm

## ICare VFit Modular Rise & Recline Lift Chair

The VFit range is a modular medical chair that offers customisable comfort with removable components, adjustable features, and innovative support technologies. It is perfect for both in-home and facility use.

### Features

- Two motors - rise and recline
- Versatile applications with removable pocket and cover
- Removable sides
- Central locking with dual wheel system allowing the chair to be moved effortlessly
- Caregiver Push Bar
- Heel Pressure Relief
- Zero gravity positioning
- Full range of cushion supports and positioning aids



### Specifications


- 170kg user weight
- Seat size - 530mm wide, 500mm deep
- Hand piece - back lit. Wired stretch cord
- Seat height - 550mm

## Updated Version of the IC333 Bed

Icare has introduced Curvex technology to its bed range, featuring rounded corners for enhanced safety and a modern design.

The curved feature protects delicate skin, shins and toes while ensuring smoother navigation and preventing wheelchair users from catching edges.


At The Mobility Store we will be stocking the IC333 bed base with the CurveX feature as well as stocking the original IC333 bed base.





## Winter Warmers

Fix the winter chills by visiting our store or website to view our range of sheepskin products designed to provide warmth and comfort for your clients. Our range includes wrap around boots, snug boots, crutch covers, protectors and overlays to stay warm this winter.



## Spacelite Ultra Light Weight Carbon Rollator

- Comfort curved backrest
- Effortless braking with locks for safety
- All carbon fibre
- User weight 150kgs
- Weighs 4.8kgs
- Comfortable ergonomic grips
- Handy utility basket
- Curb Climber

Try it  
in store  
today !



### Fun fact from TMS team

Did you know that the “winter blues” aren't just in your head? Shorter days and less sunlight can actually upset your brain’s serotonin level. The good news is getting outside for just 10–15 minutes a day can give your mood a much-needed boost.



We always love hearing from our Allied Health community – so keep the feedback coming and don't be shy about popping in for a chat or a product demonstration.

If you have any questions or would like to learn more about these products, feel free to visit our store or contact our team—we're happy to help!

Until next time,  
The Mobility Store Team

You will find The Mobility Store in 2 locations:

Unit 1, 75 Erindale Rd, Balcatta, WA ☎ (08) 6244 0033

Shop 2, 214-218 Whitehorse Rd, Blackburn, VIC ☎ (03) 9894 2405

