



**PULLS APART**  
for transport and storage

Part of the **FOR-DE Group Pty Ltd**

**NEW COLOUR!**



**PORTABLE  
SCOOTER**



METALLIC  
BLUE  
IK4SEB

METALLIC  
ORANGE  
IK4SEO

METALLIC  
GREEN  
IK4SEG

**MOVE**  
PORTABLE MOBILITY SCOOTER

**FEATURES**

- Portable scooter
- Pulls apart for transport and storage
- Comfortable mid back contoured seat
- Full suspension
- Padded loop handlebar
- Angle adjustable tiller
- Width adjustable flip up padded armrests
- Removable shopping basket
- LED Headlight
- Horn
- Fingertip operated controls
- Forward / Reverse drive settings
- Variable speed setting
- Battery power indicator
- Rear wheel drive
- Flat free tyres
- Electromagnetic brake
- Front nudge bar
- Anti tip wheels at rear

**WARRANTY**

- Structural Frame: **5 years**
- Drive Train Assembly: **18 months**
- Other Parts: **12 months**
- XT 20 Amp Hybrid Gel Batteries: **18 months**
- 2 Amp Charger: **18 months**

**Flip over for Specifications**

P: **08 9248 4180**  
E: **orders@for-de.com.au**

**ASK FOR THE BRAND**

**ikonmobility.com.au**

IMPORTANT NOTE: It is the responsibility of the user to ensure the suitability of this product for its intended use. The user should consult their appropriate healthcare professional if unsure of the correct product usage for their unique circumstances.


IKON MOBILITY strives to provide the highest quality products possible that are subject to change without notice.



Please contact IKON Mobility if you would like any further information on the Move Portable Mobility Scooter.


# MOVE

PORTABLE MOBILITY SCOOTER




OVERALL LENGTH

1080mm




OVERALL WIDTH

590mm



OVERALL HEIGHT

880mm




FRONT ASSEMBLY

L: 940mm  
W: 450mm  
H: 410mm  
Weight: 19kg




REAR WHEEL ASSEMBLY

L: 320mm  
W: 49mm  
H: 250mm  
Weight: 13kg




BATTERY PACK

W: 330mm  
D: 230mm  
H: 230mm  
Weight: 16kg




SEAT (FOLDED DOWN)

W: 460mm  
D: 610mm  
H: 530mm  
Weight: 11kg




BASKET

W: 270mm  
D: 210mm  
H: 240mm  
Weight: 1kg




SEAT WIDTH (INT)

450mm




SEAT DEPTH

440mm




SEAT HEIGHTS

330mm  
335mm  
380mm




ARMREST LENGTH

320mm




ARMREST WIDTH

80mm




ARMREST HEIGHT

200mm




TILLER WIDTH

400mm



TILLER HEIGHT

760mm



TURNING CIRCLE

1400mm



SPEED

6.4km / h



MAXIMUM RANGE

Up to 18km




SAFE CLIMB ANGLE

6°




MAXIMUM CLIMB ANGLE

10°




MOTOR

270w high torque motor



BATTERY PACK

2x 24V / 20 Ah  
(supplied)



BATTERY CHARGER

2 Amp  
(supplied)



TYRES (FLAT FREE)

Ø200mm



PRODUCT WEIGHT

60kg



MAXIMUM USER WEIGHT

130kg  
(on flat ground)

ASK FOR THE BRAND

ikonmobility.com.au

**IMPORTANT NOTE:** It is the responsibility of the user to ensure the suitability of this product for its intended use. The user should consult their appropriate healthcare professional if unsure of the correct product usage for their unique circumstances.

IKON MOBILITY strives to provide the highest quality products possible that are subject to change without notice.