

# Ultra Ryder Carbon Fibre Power Wheelchair

Lightweight and  
highly manoeuvrable



**Excellent manoeuvrability**  
makes this power  
wheelchair ideal for  
indoor environments and  
navigating local shops.



**Lightweight carbon fibre  
and aluminium frame**  
combines strength with  
durability offering a sleek,  
modern design. With a  
minimum range of 20km,  
you can go further on  
every ride.



## Effortless Mobility, Maximum Comfort

This versatile power wheelchair offers durability and ease of use. Excellent manoeuvrability makes it perfect for indoor environments and local shops, providing smooth, effortless control for everyday activities. With a lightweight and compact front folding frame the Ultra Ryder is easy to transport to family outings or holidays.



**Suspension** provides a softer and smoother ride, enhancing comfort and stability.



**Front folding frame** for effortless handling and transportation.



## Specifications

### KA390

Overall Dimension	W620xD920xH1000mm
Folded Dimension	W620xD770xH320mm
Floor to Seat Height	560mm
Seat Width	410mm
Seat Depth	420mm
Width between Armrests	465mm
Front Wheel Specification	Aluminium, Solid Tyre, 180mm
Rear Wheel Specification	Magnesium, Solid Tyre, 275mm
Net Weight (excl. battery and accessories)	17kg
Battery Weight	1.8kg
Maximum User Weight	150kg
Maximum Speed	6km/h
Battery Range*	Min. 20km

**Remote control** for attendant use, or to make independent transfers to or from the wheelchair easier.



**Swing back armrests** make it easier to transfer on or off the wheelchair.



\*(Minimum Potential Range is tested on a flat road with average speed. It will vary due to driving surface, operation habit, battery age and temperature etc.)

Please contact [customerservice@kcare.com.au](mailto:customerservice@kcare.com.au) for more information.  
[www.kcare.com.au](http://www.kcare.com.au) | 1300 783 783