



A better wheelchair.  
A better fit.



# Karma Eagle

CODE KM-8020Q

✓ DVA APPROVED

✓ TAC APPROVED

✓ AUSTRALIAN STANDARDS



<b>Frame</b>	Light weight aluminium Folding
<b>Footrest</b>	Swing away removable Height adjustable
<b>Armrest</b>	Swing back   Height adjustable Ergonomic pads
<b>Wheels</b>	Quick release   24" x 1 3/8" spoke   Pneumatic tyres
<b>Castors</b>	Solid 8"
<b>Brakes</b>	Pull to lock
<b>Seat</b>	Fire retardant upholstery Height adjustable
<b>Back</b>	Removable padding Fire retardant upholstery Height adjustable
<b>Warranty</b>	10 Years - Frame

The Karma Eagle is a lightweight, durable and reliable manual wheelchair that makes it the right choice for both users and carers.

The Eagle offers an impressive range of standard features with an extensive selection of personal customising available. The Eagle wheelchair is



easy to fold and transport and is the logical choice for people seeking a quality mobility product at an affordable price.

The use of aircraft-grade materials and unique design of the cross-bars and link bars provide this model with outstanding strength and durability.

## OPTIONS

See pages 10-11

## SPECIFICATIONS

Description	Seat Width	Seat Depth	Seat Height	Back Height	Max User Weight	Chair Weight
Eagle 16"x16"	410mm	410mm	500mm	400mm	130kg	14.4kg
Eagle 18"x16"	460mm	410mm	500mm	400mm	130kg	14.6kg
Eagle 18"x18"	460mm	460mm	500mm	400mm	130kg	14.6kg
Eagle HD 20"x18"	510mm	460mm	520mm	400mm	160kg	15.7kg
Eagle HD 22"x18"	550mm	460mm	520mm	400mm	160kg	15.9kg