

Hello from The Mobility Store!

It's hard to believe we're already into April. Time certainly flies, and we want to take a moment to express our sincere gratitude for the continued support from our wonderful community of Allied Health Professionals. You make all the difference!

What's New at The Mobility Store?

As we look ahead in 2025, we're thrilled to share that we are expanding our range of scooters and powered wheelchairs. We're always striving to provide the best mobility solutions, and these exciting new additions are sure to make a difference in the lives of your clients.



Scooter Stopper

The Scooter Stopper is an ideal accessory to give a support person/carer for safe stopping control over a scooter/wheelchair user in any situation.

Utilises a modern, reliable remote control to safely bring the user to a stop from a distance of 25 to 45 metres, depending on the working environment.

This device is powered directly through the equipment's charger socket and ensures the user cannot operate the machine until the remote is reset. Tested on most popular scooters and powered wheelchairs, it's a fantastic tool for added safety and peace of mind.

Explore the New Pride Range of Scooters

We're also excited to showcase a **NEW** selection of Pride Mobility Scooters in our showroom! We recently introduced the Go Go Carbon Scooter, the Baja Scooter (with a 185kg weight capacity), and the 130XL to our range.

Pride Mobility is known for their exceptional quality.

Go Go Carbon Fibre Features

- 136 kg weight capacity
- 134.6 cm turning radius
- Up to 6 km/h maximum speed
- Up to 16 km per charge
- Folds down in just seconds and can be towed like luggage
- Full carbon fibre frame, rims and tiller
- NFC ignition provides extra security
- LED lighting for enhanced visibility
- Removable back seat and base for easy cleaning
- Include turning speed sensor to prevent turns at excessive speeds



Pathrider 130XL Features

- 158.7 kg weight capacity
- 167.6 cm turning radius
- Up to 10 km/h maximum speed
- Up to 24.9 km per charge
- Large 33cm pneumatic tyres for excellent outdoor performance
- Full lighting package including directional signals and angle-adjustable, dual bulb lower LED headlight
- Seat post suspension for a smoother ride
- Easy drive tiller with wraparound handles
- Easily accessible tie down points (for transport of unoccupied scooter)



Baja Bandit Features

- 181 kg weight capacity
- 134.7 cm turning radius
- Up to 10 km/h maximum speed
- Up to 23.3 km per charge
- Front, rear & seat post suspension
- Angle adjustable tiller
- Front & rear bumper guard



Robooter Powered Wheelchair Collection

We're excited to add the stylish and innovative Robooter power wheelchairs to our collection, offering seamless transitions between outdoor adventures and indoor comfort.



Robooter E40 Features

- Folding mode manual
- Maximum speed 8km
- Maximum barrier height 4cm
- Maximum crossing width 10cm
- Slope climbing performance 9 degrees
- Total weight 25kg
- Battery weight 3.4kg
- Maximum load 150kg
- Battery capacity 20ah/25ah
- Overall dimensions 100x62.4x93cm
- Folding shape dimensions 70x62.4x45cm
- Seat width 48cm
- Travel distance 23km-29km depending on battery size

Robooter E60 Features

- Folding mode manual
- Maximum speed 8km
- Maximum barrier height 6cm
- Maximum crossing width 10cm
- Slope climbing performance 15 degrees
- Total weight 32.4kg
- Backrest angle, not adjustable 18 dgs
- Battery weight 3.3kg
- Specifications of front and rear wheels 10 inch
- Maximum load 150kg
- Battery capacity 20ah/25ah
- Overall dimensions 98.7x64x93.2cm
- Folding shape dimensions 77x64x50cm
- Seat width 48cm
- Travel distance 23km-29km depending on battery size



Thank you for your continued support.

We're excited about these new additions and, as always, grateful for your ongoing trust in us.

If you have any questions or would like to learn more about these products, feel free to visit our store or contact our team—we're happy to help!

Until next time,
The Mobility Store Team

Stay connected with The Mobility Store

