








Ellipse Lite Aluminium Rollator

Product Code: 9346376000536 (Small, Red)

Features:

- AU STD. Tested 492485
- Adjust Height Tool-Free
- Adjust Brake Tool-Free
- Foldable Waterproof Seat
- Easily Installable Nylon Bag



| | | | | | |
|---------------------|--|----------------------------------|--|----------------------------|---|
| Length (cm) |  59.5 | Folded Width (cm) |  26 | Front Castor Diameter (cm) |  24.2/9.5" |
| Width (cm) |  49 | Seat Depth (cm) |  22 | Rear Castor Diameter (cm) |  20.3/8" |
| Max Height (cm) |  82 | Seat Width (cm) |  45.4 | Max Load (kg) |  136 |
| Minimum Height (cm) |  73.5 | Seat Height from the Ground (cm) |  49 | Net Weight (kg) |  8 |