Unit 11/613 Whitehorse Road
Mitcham
VIC 3132
Australia


Ph: (03) 9874 5785
www.quintro.com.au



ALPHA 428 ROLLATOR

Features:

- lightweight aluminium frame
- comfortable 2" padded seat
- easily removable foam backrest
- adjustable handle & loop brakes
- 8" wheels & nylon storage bag
- easy to fold for travel or storage
- user weight capacity of 125kg



Alpha 428 Rollator	**Metric**	**Imperial**
Approximate Total Product / Unit Weight	8.5 kg	18.7 lbs
Overall Width (product or frame when opened)	630 mm	24.8 inch
Overall Depth (product or frame when opened)	700 mm	27.6 inch
Width / Depth (product or frame when closed)	510 mm	20.1 inch
Overall Height (or height range if adjustable)	835 - 970 mm	32.9 - 38.2 inch
Armrest Width (outside dimension)	N / A	N / A
Armrest Width (inside dimension)	440 mm	17.3 inch
Armrest Height (from top of seat)	N / A	N / A
Back / Back-Rest Height (from top of seat)	225 mm	8.9 inch
Seat Height (from top of seat to floor)	570 mm	22.4 inch
Seat Width (actual seat, not available space)	350 mm	13.8 inch
Seat Depth (actual seat, not available space)	320 mm	12.6 inch
Wheels (dimensions of front)	203 mm	8.0 inch
Wheels (dimensions of rear)	203 mm	8.0 inch
Recommended Maximum User Weight	125 kg	275.6 lbs

