ShowerGlyde
MODEL SG1










Model SG1 Nothing Adjacent to Tub



ShowerGlyde Model SG1 is a revolutionary sliding transfer bench. It is designed to fit in environments where there are no obstructions adjacent to the tub (i.e. vanity, toilet). SG1 slides easily over the tub and locks in two points; on the inside and outside of the tub. This allows for safe transfers while using the bench. SG1 can be assembled as a left or right-hand unit. One grab bar installed at the wall side of the rails is included as an aid for sliding. The user can pull or push themselves easily along the rails. SG1 is height adjustable permitting clearance of up to a 49.5cm tub height. Optional extension legs are available for taller tubs (up to 63cm high).

Model SG1 also includes a Seat Belt, a Grab Bar, and a Commode Seat with 3 options; an open Commode Seat for easy washing, a Reducer Ring to reduce the size of the commode opening, and a Solid Seat Insert to close the commode opening completely.

SWL 150KG

Includes:


- Grab Bar
- Seat Belt
- 49.5 cm tub clearance
- Commode Seat
- Reducer Ring Insert
- Sold Seat Insert


showerbuddy®


ShowerGlyde


MODEL SG1

SolutionBased










Tub clearance.

Add additional clearance
with leg extension
accessory: PD000-298
502-640 mm 19.75"~25"



For more information contact
admin@shower-buddy.com
www.shower-buddy.com