# ShowerGlyde

MODEL SG2










## Model SG2
## Toilet Adjacent




ShowerGlyde Model SG2 is a sliding transfer bench designed to fit in environments with a toilet adjacent to the tub. Similar to SG1, SG2 slides effortlessly over the tub and locks in two points; on the inside of the tub and over the toilet. A longer set of rails accommodates the use of the commode over the abutting toilet. SG2 can be assembled as a left or right-hand unit. One grab bar installed at the wall side of the rails is included for support in sliding. SG2 is height adjustable allowing clearance of up to a 49.5cm tub height.
Optional extension legs are available for taller tubs (up to 63cm high).

Model SG2 also includes a Seat Belt, Grab Bar, and a Commode Seat with 3 options; an open Commode Seat for over-toilet use, a Reducer Ring Insert to reduce the size of the commode opening, and a Solid Seat Insert to close the commode opening completely.

SWL 150KG

### Includes:


- Grab Bar
- Seat Belt
- 49.5 cm tub clearance
- Commode Seat
- Reducer Ring Insert
- Sold Seat Insert


showerbuddy®


# ShowerGlyde


MODEL SG2

SolutionBased





Add additional clearance with leg extension accessory: PD000-298
502-640 mm 19.75"~25"




For more information contact
admin@shower-buddy.com
www.shower-buddy.com